- Compensation for a metabolic acidosis is hyperventilation to decrease the arterial pCO2.
- The metabolic acidosis is detected by both the peripheral and central chemoreceptors and the respiratory center is stimulated. The initial stimulation of the central chemoreceptors is due to small increases in brain ISF [H+]. The subsequent increase in ventilation causes a fall in arterial pCO2 which inhibits the ventilatory response.
 - Maximal compensation takes 12 to 24 hours
 - The chemoreceptor inhibition acts to limit and delay the full ventilatory response until bicarbonate shifts have stabilised across the blood brain barrier. The increase in ventilation usually starts within minutes and is usually well advanced at 2 hours of onset but maximal compensation may take 12 to 24 hours to develop.
- The arterial pCO2 at maximal compensation has been measured in many patients with a metabolic acidosis. A consistent relationship between bicarbonate level and pCO2 has been found. It can be estimated from the following equation:

Expected pCO2 = 1.5 (Actual [HCO3]) + 8 mmHg (Units: mmols/I for [HCO3], and mmHg for pCO2).

- The limiting value of compensation is the lowest level to which the pCO2 can fall this is typically 8 to 10mmHg, though lower values are occasionally seen.
- If a patient with a severe metabolic acidosis requires intubation and controlled ventilation in hospital, the acidosis can markedly worsen unless the hyperventilation is maintained. The ventilation should be set to mimic the compensatory hyperventilation to keep the pCO2 low.
- Carbon dioxide crosses cell membranes readily so intracellular pH falls rapidly also, resulting in depression of myocardial contractility, arrhythmias and a rise in intracranial pressure.

General principles are as follows:

- Accurate diagnosis of the cause of the metabolic acidosis is essential because this allows correct treatment of the underlying disorder
- 2. Treat the underlying disorder as the primary therapeutic goal
- Provide supportive treatment (eg fluids, oxygen, treatment for hyperkalaemia) including all appropriate emergency management
- 4. In most cases, IV sodium bicarbonate is NOT necessary, NOT helpful, & may even be harmful in the treatment of metabolic acidosis.

Kidney: Renal generation of new bicarbonate

- This usually occurs as a consequence of an increase in ammonium excretion.

Liver: Hepatic metabolism of acid anions to produce bicarbonate

- The normal liver has a large capacity to metabolise many organic acid anions (eg lactate, ketoanions) with the result that bicarbonate is regenerated in the liver. In severe ketoacidosis there is often a large loss of ketoanions due to the hyperglycaemia induced osmotic diuresis. This leaves a shortfall of ketoanions to be used to regenerate bicarbonate as a consequence of their metabolism in the kidney.

Exogenous Administration of sodium bicarbonate

- This is the time honoured method to 'speed up' the return of bicarbonate levels to normal. Indeed, this may be useful in mineral acidosis (hyperchloraemic metabolic acidosis) where there are no endogenous acid anions which can be metabolised by the liver. However, in most other cases of metabolic acidosis this administration is either not helpful or may be disadvantageous.
- A metabolic acidosis is often strongly suspected because of the clinical presentation of the patient (eg diabetes, renal failure, severe diarrhoea). Three clues from a typical hospital automated biochemical profile are:
- (i) Low 'bicarbonate' (or low 'total CO2')
- (ii) High chloride
- (iii) High anion gap
 - In addition to arterial blood gases, some other investigations useful for indicating a metabolic acidosis and for differentiating between the various major causes are:
 - (i) Urine tests for glucose and ketones
 - (ii) Electrolytes (incl chloride, anion gap, 'bicarbonate')
 - (iii) Plasma glucose
 - (iv)Urea and creatinine
 - (v) Lactate
- useful additional indices in assessment of metabolic acidosis include:
- (i) Anion gap
- (ii) Delta ratio
- (iii) Urinary anion gap
- (iv) Osmolar gap

